# TISSIDental

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INOX - WHITE - PLUS - 4 BAR - 6 BAR

**Manuale di Istruzioni**  
**Instruction Manual**  
**Gebrauchs Anweisungen**  
**Mode d'emploi**  
**Manual de manejo**

**CAUTION**

INSTALL THIS DEVICE IN A PROPER LOCATION  
TO BE USE BY AUTHORIZED AND INSTRUCTED PERSONNEL ONLY

# INDICE / INDEX

Declaration of conformity	3
Exploded draw	4
Spare parts list	4
Electrical scheme	6
Visual instructions	7
Water CAUTION	7
Safecap	8

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---

## ENGLISH

Instruction Manual	12
Safety directives	12
Unpack	12
Installation	12
Maintenance	13
Guide to common questions	13
Specifications	13
Dimensions	13

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TISSI reserves the right to change the specifications of this equipment without notice.

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Herewith we declare that the product is in compliance with the relevant requirements in the following directives:  
 2006/42/EC (Machinery safety)  
 2006/95/EC (Low voltage equipment)  
 2004/108/EC (Electromagnetic compatibility)

Harmonized specifications applied:  
 EN ISO 12100-1:2010, EN 61010-1:2010,  
 EN 61326-1:2006

Authorised to compile the technical documentation:  
 Tissi G. I.

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

- PL WIRE PLUG
- PD FOOT PEDAL
- EV ELECTRICAL STEAM VALVE
- LH LAMP HOLDER WATER LEVEL
- HL LAMP HEATING ELEMENT
- PS PEDAL PANEL SWITCH
- GS GENERAL SWITCH
- TH THERMOSTATIC CONTROL
- PR PRESSURE CONTROL
- HE HEATING ELEMENT 1,5 KW
- CN CONDENSATOR
- GR GROUND
- FS FUSE 10 AMP – T/10A 500V 6,3x32



**ELECTRICAL SCHEME PLUS MODEL (GUN WITH SWITCH)**  
**SCHEMA ELETTRICO MOD PLUS (PISTOLA CON INTERRUETTORE)**

- GR ground – terra
- FS fuse 6x32 20 a/t
- GS gen switch – int generale
- CN radio filter – filtro radio
- GL gen lamp - spia generale
- HE heating coil – resistenza
- R rele' heat – rele' resistenza
- RVS e.valve switch – int e. valvola
- PVL int light switch – int luce int
- RV rele' e.valve – rele' e. valvola
- LV internal light – luce interna
- T transformer – trasformatore
- FV fuse 6x32 1 a/t
- EV electrical valve – elettrovalvola
- MS microswitch – microinterruttore

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
## VISUAL INSTRUCTIONS

*Please read instruction booklet carefully before installing and use of this device*

Read instruction booklet before operate

Connect steam gun hose

Insert gun holder

Insert plug (check voltage)

Check NO pressure


Switch ON left button

Slowly, unscrew the cap

(Press the switch on the gun with the PLUS model)

Insert funnel and fill with demineralized water to the MAX level

Switch OFF left button

Switch ON right button (MAIN)

## WATER CAUTION

- **ATTENTION: USE ONLY WATER WITH LOW MINERALS CONTENT (DISTILLED OR DEMINERALIZED WATER)**
- **ATTENZIONE: UTILIZZARE SOLO ACQUA A BASSO CONTENUTO DI SALI MINERALI (ACQUA DISTILLATA O DEMINERALIZZATA)**
- **AUFMERKSAMKEIT: VERWENDEN SIE NUR DESTILLIERTES ODER ENTMINERALISIERTES WASSER**
- **ATTENTION: EMPLOYEZ SEULEMENT L'EAU DISTILLÉE OU DÉMINÉRALISÉE**
- **ATENCIÓN: UTILICE SOLAMENTE AGUA BAJA DEL CONTENIDO MINERAL (DESTILADA O DESMINERALIZADA)**

<p><b>NO</b></p> 	<p><b>OK</b></p> 	<p><i>Non stringere il tappo con forza</i></p> <p><i>Do not overtight the stopper/safety valve</i></p> <p><i>Das sicherheitsventil nicht zu fest aufschrauben</i></p> <p><i>Ne pas serrer excessivement la soupape de sûreté</i></p> <p><i>No apriete en exceso el tapón-válvula seguridad</i></p>
--	--	--

The steam generators are equipped, in the greatest part of the cases, of simple stoppers with integrated a safety valve. This safety valve only enters in function when a certain pressure has been reached inside the boiler, releasing controlled steam towards the outside. This avoids to reach dangerous pressures inside the boiler, for its own integrity. These stoppers do not represent a device of emergency for the customer, because they do not prevent the removal when the boiler is in pressure. Therefore “SAFECAP” is the new patented stopper that is at the same time safety valve and antiremoval device. When the pressure in the container arises over 0,3 bars, it is not possible to unscrew the stopper, while when it is necessary the action of the safety valve, “SAFECAP” allows the release of the steam towards the outside in one not dangerous direction: that is towards the bottom. The internal components are in stainless steel and the gaskets are made with a new material with a very high reliability. The external gasket is in Viton and can be periodically substituted. We advise to change the SAFECAP every 18 months to increase safety. SAFECAP has been designed to perform two tasks:

- 1 to prevent the opening of the boiler through the removal of the same stopper, when in the boiler there is a pressure over 0,3 bars. When it is 0,3 bar ( $\pm 0,05$  bar), the grip handle comes unhooked from the body of the stopper and turns to empty preventing its removal.
- 2 in case of inside overpressure to allow the controlled steam with the jet towards the bottom. When the steam pressure (5,5 bars  $\pm 0,5$  bar) has been reached, the steam will be released through the openings in the lower part of the grip handle. The removal of the stopper only returns possible when the pressure inside the boiler comes down newly under 0,3 bar.

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## CONGRATULATIONS FOR YOUR CHOICE!!

Choosing this apparatus you have purchased a particularly reliable unit that will facilitate your work.

### IMPORTANT !

**Before installing the apparatus read carefully this manual and keep it for future consultations.**

### WARNING

Do follow carefully all the indications brought back in this manual. During all the maintenance operations disconnect the power cord from the electrical outlet, and/or disconnect the compressed air and/or the water attack and/or other attacks of external parts to the apparatus. Do not introduce stranger objects inside the unit during the maintenance operations. Use only the accessories in equipment for the filling operations and/or maintenance. Do not execute different operations of maintenance from those brought back in the manual. Whichever operation not contemplated in this manual, can involve risks. For operations of this type, contact the technical service. Before turning on the unit and connecting the power cord to the electrical outlet, make sure that the voltage corresponds to the one brought back on the fabrication nameplate placed on the back part of the apparatus.

### CAUTION

- **Do not make any maintenance operation when the power cord is connected: FULGURATION DANGER!!!**
- **Use only genuine spare parts. Placing not genuine spare parts can void your warranty, moreover can make the apparatus not in compliance with the norms, and dangerous when in use.**
- **Do not touch the electrical valve body – WARNING: HOT SURFACE**
- **For whichever inherent information concerning installation, maintenance and use of the unit, contact the customers' service TISSI:**

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*Fax (+39) 02 98285327 - Tel. (+39) 02 98282694 - www.tissidental.com*

### UNPACK

Verify integrity of the pack at the moment of the reception of the apparatus. Verify on the back side the yellow adhesive nameplate QUALITY CONTROL PASSED the control code is integral. Verify that in the packaging all the parts and the supplied accessories of series are contained together with the apparatus: steam gun with a set of 4 nozzles, loading funnel water, rod rest-gun, electrical pedal, user's documentation and warranty.

### INSTALLATION

Take the unit out of the packaging and place it on a flat surface strong enough to bear its weight (14 Kg. with water). Connect the electrical pedal (ref. 2) to the appropriate outlet placed on the upper side of the unit. Insert the rod rest-gun and screw it in the appropriate site (ref. 40). Connect the steam tube (ref. 1) to the appropriate sleeve of electric valve (ref. 9) and properly fix it with the equipped band. Connect the main plug to a power outlet.

**ATTENTION:** make sure that the voltage corresponds to the one brought back on the fabrication nameplate placed on the back side of the apparatus. Press the power switch on the left to position 1 (ref. 31 SX) to turn on the unit and the inside water level light. Remove the filling stopper (ref. 8) making a light pressure downwards; then fill the boiler with demineralized water, until you reach the maximum level showed on the front panel of the unit. Screw the stopper (ref. 8) on the hole of threaded loading to close. Now turn off the left switch and turn on the switch on the right (ref. 31 DX) to start the heating; the ignition of the heating element is marked by the red light (ref. 33) placed on the left side of both the two switches. When the pressure of 3.5 bars is reached the red

light automatically switch off, now the unit is ready to work. Possible boiling noise when starting can be eliminated letting out the inside air bubbles by pressing the steam spilling pedal. To obtain the steam out of the gun (ref. 1), turn on the switch on the left (the same which turns on the water level light), then press the electrical pedal (ref. 2) previously connected. Regulate the steam flow by rotating the regulation grip handle placed on the base of the electric valve (ref. 9). The initial steam capacity could be very short, wait for a few minutes till the pressure is reached again, then the steam will be released continuously. When you end steaming, turn off the switch on the left to disconnect the pedal and electric valve. At the end of the day turn off the main switch on the right, without unloading the pressure.

## MAINTENANCE

Check the level of the water: verify every time the unit is turned on that the level of the water is not under the minimum level reference showed on the front panel. If the level of the water is too low, fill the boiler following the procedure of loading of the water (see points 6-7-8).

- **Drainage of the residual calcium in the boiler** - Every 6 months of operating, the drainage of the residual calcium in the boiler must be carried out executing the following operations: control that the apparatus is cold and without internal pressure. Disconnect the power cord from the electrical outlet (rif. 39). Place the apparatus on a bench on its left side. Unscrew the nut placed under the base of the boiler (rif. 54) without removing it completely to avoid to bathe the inside of the apparatus. Place under the apparatus a washbasin or a container of at least 5 liters capacity. Then unscrew completely the nut from the boiler and remove it from the center. Remove now the stopper (rif. 8) and make discharge the residual water completely. Pour inside with the funnel 3/4 liters of water for rinsing the boiler. Then screw again not too strongly the nut under the boiler, placing a tape of teflon on the thread (approximately 8/9 complete turns). Reconnect now the apparatus to the power cord, and fill the boiler with demineralised water as usual.

## GUIDE TO COMMON QUESTIONS

### The apparatus is not turning on:

- 1 Control that the power cord plug is properly connected to the electrical outlet.
- 2 Check that power is on at the electrical outlet.

### Steam flows little in relation to the pressure indicated on the gauge:

Check that the regulation grip handle of the steam placed on the electric valve (rif. 9) is completely opened.

### A strong water boiling noise from the boiler when start heating (rare operation):

Open slowly the filling cap to let air bubbles out.

### After reached the pressure the steam comes out for limited time:

The machine's components are cold, please allow necessary time to warm up all internal components (15-20 minutes at 3.5 bar without release of steam) Carried out these controls, in case the apparatus still presents disadvantages, contact the customers' technical service.

## SPECIFICATIONS

W/Hz 230/50 - 60 or 110/50  
Absorption 1500 Watt (1300W 110V)  
Weight Kg. 10

## DIMENSIONS

Height cm. 50 - Width cm. 22 - Depth cm. 30

## NOTES

- For the external cleaning of the device use a cotton tawel damped with water.
- The device can work in laboratory room from 0°C to 35°C.
- Noise level (1 mt): stand-by 26.6 db, warm-up 37.5 db, steam cleaning 68.7 db.
- Maximum Pressure in Bar/Hpa (Kilopascal): 3,5/350

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