TISSIDental

WATER JET

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CAUTION

INSTALL THIS DEVICE IN A PROPER LOCATION
TO BE USE BY AUTHORIZED AND INSTRUCTED PERSONNEL ONLY

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ENGLISH

Instruction Manual	8
Warning	8
Unpack	8
Installation	8
Technical data	10
Maintenance	10

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TISSI reserves the right to change the specifications of this equipment without notice.

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Herewith we declare that the product is in compliance with the relevant requirements in the following directives:

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CONGRATULATIONS FOR YOUR CHOICE!!

Choosing this apparatus you have purchased a particularly reliable unit that will facilitate your work.

IMPORTANT !

Before installing the apparatus read carefully this manual and keep it for future consultations.

WARNING

- Do follow carefully all the indications brought back in this manual.
- During all the maintenance operations disconnect the power cord from the electrical outlet, and/or disconnect the compressed air and/or the water attack and/or other attacks of external parts to the apparatus.
- Do not introduce stranger objects inside the unit during the maintenance operations.
- Use only the accessories in equipment for the filling operations and/or maintenance. Do not execute different operations of maintenance from those brought back in the manual.
- Whichever operation not contemplated in this manual, can involve risks. For operations of this type, contact the technical service.
- Before turning on the unit and connecting the power cord to the electrical outlet, make sure that the voltage corresponds to the one brought back on the fabrication nameplate placed on the back part of the apparatus.

CAUTION:

- **Do not make any maintenance operation when the power cord is connected: FULGURATION DANGER!!!**
- **Use only genuine spare parts. Placing not genuine spare parts can void your warranty, moreover can make the apparatus not in compliance with the norms, and dangerous when in use.**
- **Do not touch the electrical valve body – WARNING: HOT SURFACE**
- **For whichever inherent information concerning installation, maintenance and use of the unit, contact the customers' service TISSIDental:**

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Fax (+39) 02 98285327 - Tel. (+39) 02 98282694 - www.tissidental.com

UNPACK

- Verify integrity of the pack at the moment of the reception of the apparatus.
- Verify on the back side the yellow adhesive nameplate **QUALITY CONTROL PASSED** the control code is integral.
- Verify that in the packaging all the parts and the supplied accessories of series are contained together with the apparatus: steam gun with a set of 4 nozzles, loading funnel water, rod rest-gun, electrical pedal, user's documentation and warranty.

USE

Device to produce steam, is ideal for cleaning operations during the construction of prosthesis, to clean fixed and mobile prosthesis after polishing, to eliminate wax from models and flasks, as well as during specific operative processes required by certain materials used for dental purposes.

INSTALLATION

Make sure the device does not show evident signs of damage due to transport.

Place the device on an horizontal surface or hang it to the wall using the hooks located on the back of the device, and the two plugs supplied with it. These plugs should be fixed to the wall 28 cm apart. **The device must not be installed above sources of heat (i.e. heaters) or above sinks.** Also make sure that the openings on the lid are not covered thus obstructing ventilation.

The steam generator can be connected directly with the water system or with an outer tank.

■ **Water system connection**

If the apparatus must be connected with the water system, it will be necessary to supply the steam generator with a water softener in order to avoid limestone sediments. A water softener is included with the machine (check it periodically and substitute after 3 years of use).

The connection with the water softener must be performed by skilled and authorised personnel, who must previously carry on the setting of the water softener in accordance with the instructions enclosed to the same softener. The non-observance of these instructions could seriously compromise the steam generator.

Connect the steam generator with the water softener by means of the pipe (with smaller diameter) supplied together with the apparatus, by introducing it into the water inlet fitting placed on the right side of the steam generator.

Connect the water softener with the water system tap by means of the pipe supplied (with bigger diameter).

Open the water system tap.

■ **Outer tank connection**

- If the steam generator must be connected with an outer tank (it is advisable to use a container with a capacity of 5 litres at least), introduce one end of the pipe (with smaller diameter) supplied together with the apparatus into the water inlet fitting placed on the right side of the steam generator and the other end into the tank previously filled with distilled water, adding a teaspoonful of table salt. To avoid overheating the pump, place the container at a minimum height of 50 cm with respect to the steam generator.

Before installing the device, make sure that the power of the supply mains complies with the rating plate's indications (230 Vac, 50/60 Hz.) and that the line is protected by a magnetothermal differential switch. Also verify the absorption and the installed power. In case of doubt about the type of available power, please contact your Electric Power Company.

Make sure that the main switch is on the (0) off position, then connect the plug of the power cable to an outlet. Make sure that the power system has a grounding circuit.

OPERATION

To use the steam generator, please carefully follow the instructions reported here below and in the proper sequence:

- 1 **Turn on the device**
- 2 **Activate pump/ resistor**
- 3 **Steam/water output**

1 TURNING ON THE DEVICE

Turn on the device by bringing the main switch on the ON (I) position.

The green light inside the switch lights up.

2 ACTIVATE PUMP-HEATING

Activate the pump by bringing the resistor/pump switch on the ON (I) position.

The green light inside the switch lights up

If the steam generator is **connected with an outer tank**, open the valve placed on the right side of the generator and always leave it open. When the pump starts working, a small quantity of water flows out of the valve. Therefore, it is advisable to place a small container under the valve. The low water level light lights up for a few seconds. When the light goes off, the resistor active light turns on. The resistor remains active until a pressure of 6 bar is reached. The operator may anyway use the device even if the pressure of 6 bar has still not been reached checking the boiler's pressure on the gauge. The operator may check if the heating is working by seeing if the light below the heating switch is lit.

3 STEAM-WATER OUTPUT

To make the steam come out, press the button on the left of the gun.

To make the water come out, press the button on the right of the gun.

Be careful not to aim the steam coming out of the gun directly towards people, animals, things and /or sources of electricity or the device itself. The gun should be used at such a distance from the device so that the steam jet does not penetrate the openings on the top and on the bottom of the steam generator.

The intensity of the steam jet may be adjusted with the knob located in front of the device:

- **clockwise rotation = flow decrease**
- **counter-clockwise rotation = flow increase**

TECHNICAL DATA

Dimensions:	W 46 x D 32 x H 35,5 cm
Weight:	Kg 15
Voltage:	230 Vac 50/60 Hz.
Power:	2200 Watts
Fuses:	6x32 mm F 10 A/250 V (fuses power supply) 5x20 mm F300 mA/250 V (fuse transformer)
Use:	Indoors
Temperature:	from +5°C to +40°C
Relative Humidity:	80% maximum
Altitude:	up to 2000 meters
Installation category:	II
Degree of pollution:	2

REGULAR MAINTENANCE

The steam generator does not require any particular maintenance. It is important to use always a water calcium cleaner (demineralized water). A device for this purpose is included with the machine.

WARNINGS

- The steam nozzle is already installed on the gun when the steam generator is delivered. If the operator should need to replace the nozzle, such operation must only be carried out with the device turned off, the main switch on OFF (0) and the power cable disconnected. Also make sure that the black gasket of the steam output brass pipe is correctly fitted and that the nozzle, when fitted, is in the locked position.
- **In case the safety valve should come open, immediately turn off the resistor switch, bleed the boiler pressure through the gun, bring the main switch on OFF (0), disconnect the power cable and call the manufacturer.**
- **If the steam generator is connected with the water system, make sure that the water pressure does not exceed 4 Bar. In this case, connect a pressure reducer calibrated to 2 bar (max.)**
- **If the steam generator is connected with an outer tank, always make sure that the pipe end is immersed into the water and that there is always water in the container.**
- To clean the device, use a dry cloth and liquid detergent (do not use flammable liquids) then dry with a soft rag. **Do not put the device, for any reason whatsoever, under running water.** Cleaning operations must be carried out while the main switch is on the off (0) position and the power cable is disconnected. Before using any type of cleaning method, other than the one recommended by the manufacturer, the user must contact the manufacturer to make sure that the method he wants to use is not likely to damage the device.
- Do not put any object on the power cables and do not place the cable in areas where it could get trampled upon.
- Disconnect the power cable if the device should not be used for a while.
- Use the original packaging if you should need to transport the device.

INSTRUCTIONS FOR DISPOSAL

This product shall not be disposed of as a mixed municipal rubbish. For the disposal comply with the provisions of the rules in force.

IMPORTANT INFORMATION

- The Manufacturing Company declines any responsibility for any apparatus installed, utilized, modified, repaired and/or altered by non competent and unauthorized personnel.

- The apparatus has been designed and built only for the uses contemplated, consequently the Manufacturing Company refuses any responsibility for damages to persons, animals, things, and the machine itself caused by improper use of this apparatus.
The improper use of this machine is dangerous!!
- Subject the apparatus to a complete annual revision by the Manufacturer.
- To clarify any doubt and/or technical explanation please contact our Technical Department.

GUIDE TO COMMON QUESTIONS

The apparatus is not turning on:

- Control that the power cord plug is properly connected to the electrical outlet.
- Check that power is on at the electrical outlet.

NOTES

- For the external cleaning of the device use a cotton tawel damped with water.
- The device can work in laboratory room from 0°C to 35°C.
- Noise level (1 mt): stand-by 26.6 db, warm-up 37.5 db, steam cleaning 68.7 db.
- Maximum Pressure in Bar: 5,5

www.imq.it

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CISQ is a member of

THE INTERNATIONAL CERTIFICATION NETWORK
www.iqnet-certification.com

*IQNet, the association of the world's first class
certification bodies, is the largest provider of management
System Certification in the world.*

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E' CONFORME ALLA NORMA / IS IN COMPLIANCE WITH THE STANDARD

ISO 13485:2016

TISSIDental

Dental, Goldsmith and Laboratory Equipment Made in Italy

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